

# Command Index

## Input

- [Processed NMR data](#)
  - [NMRPipe](#)
  - [Felix](#)
  - [Varian VNMR](#)
  - [Bruker XWinNMR and UXNMR](#)
- [UCSF NMR format](#)
- [Sparky spectrum annotation \(aka "save"\) files](#)
- [Sparky project files](#)

## Output

- [Peak lists](#)
- [Chemical shift list](#)
- [Chemical shift table](#)

## File commands

<a href="#">file open...</a>	fo open UCSF NMR data, a Felix matrix, or a Sparky spectrum file
<a href="#">file save</a>	fs save a Sparky spectrum annotation file
<a href="#">file save as...</a>	fa save spectrum with a new name
<a href="#">project open...</a>	jo open several spectra at once
<a href="#">project save</a>	js save all spectra
<a href="#">project save as...</a>	ja save project with a new name
preferences...	pf miscellaneous settings
quit	qt

## View commands

<a href="#">axes in Hertz</a>	xh display view scale axes in Hz
<a href="#">axes in ppm</a>	xp display view scale axes in parts per million
<a href="#">axes show index</a>	xd show data index values on view scales
<a href="#">axis nucleus type</a>	xa show nucleus type (1H, 13C, 15N) on axes
<a href="#">axis roll</a>	xr exchange axes keeping center point of fixed
<a href="#">axis transpose</a>	xx zoom to cross-diagonal region of spectrum
<a href="#">center view on peak</a>	vc center a view on the selected peak
<a href="#">orthogonal views</a>	ov show 3 orthogonal views of a 3-D spectrum
<a href="#">region goto</a>	rg zoom to a named region
<a href="#">region tool...</a>	rt name regions of a spectrum

<a href="#">renumber view</a>	rv renumber view name to use lowest available suffix number
<a href="#">show axis scales</a>	vs show scales (ppm, hz, or index units) along view axes
<a href="#">show grids</a>	gs toggle display of grids
<a href="#">show labels</a>	ls toggle display of labels
<a href="#">show lines</a>	is toggle display of lines
<a href="#">show ornaments</a>	os toggle display of all ornaments
<a href="#">show peaks</a>	ps toggle display of peaks
<a href="#">show peak groups</a>	xs toggle display of peak groups
<a href="#">show peak info</a>	vp show peak position / height / linewidth in status bar
<a href="#">show resonances</a>	vR show resonances along edge of spectrum
<a href="#">show scrollbars</a>	vb show scrollbars for view windows
<a href="#">show slices</a>	vS show 1-D data cross-sections along edge of spectrum
<a href="#">view duplicate</a>	vd make another view window
<a href="#">view hide</a>	vh undisplay a view
<a href="#">view redraw</a>	vr redraw the view (in case it gets messed up)
<a href="#">view settings...</a>	vt set crosshairs, aspect ratio, slices, ...
<a href="#">zoom full</a>	zf show entire spectrum
<a href="#">zoom in</a>	zi zoom in by a factor of 2
<a href="#">zoom out</a>	zo zoom out by a factor of 2
<a href="#">zoom down</a>	zd move down a plane in a 3D spectra
<a href="#">zoom up</a>	zu move up a plane in a 3D spectra
<a href="#">zoom previous</a>	zp zoom to the previously shown region
<a href="#">zoom next</a>	zn reverse of zoom previous (zp) command

## Peak commands

<a href="#">add sweepwidth</a>	a1 a2 a3 a4	adjust peak w1,w2,w3,w4 frequency one sweepwidth downfield
<a href="#">subtract sweepwidth</a>	A1 A2 A3 A4	adjust peak w1,w2,w3,w4 frequency one sweepwidth upfield
<a href="#">reflect upfield</a>	u1 u2 u3 u4	reflect peak w1,w2,w3,w4 frequency about upfield boundary
<a href="#">reflect upfield</a>	d1 d2 d3 d4	reflect peak w1,w2,w3,w4 frequency about downfield boundary
<a href="#">assignment copy</a>	ac	remember peak assignments for assignment paste
<a href="#">assignment delete</a>	aD	delete the assignments of selected peaks
<a href="#">assignment paste</a>	aP	paste assignments onto peaks in this spectrum
<a href="#">assignment paste &amp; label</a>	ap	copy assignments and show labels
<a href="#">crossdiagonal assignment copy</a>	aX	copy assignments from crossdiagonal peaks
<a href="#">assignment tool...</a>	at	make assignments
<a href="#">delete unused resonances</a>	dr	delete resonances not used in any peak assignments
<a href="#">integration tool...</a>	it	select integration method
<a href="#">sum over box method</a>	ib	switch to box integration method

<a href="#">sum over ellipse method</a>	ie	switch to ellipse integration method
<a href="#">Gaussian fit method</a>	ig	switch to Gaussian fit integration
<a href="#">Lorentzian fit method</a>	il	switch to Lorentzian fit integration
<a href="#">integrate selected peaks</a>	pi	
<a href="#">ornament copy</a>	oc	record ornaments for future paste operation
<a href="#">ornament paste</a>	op	paste ornaments into this spectrum
<a href="#">crossdiagonal ornament copy</a>	oX	copy ornaments across diagonal of 2D homonuclear spectrum
<a href="#">center peaks</a>	pc	center selected peaks at local maxima or minima.
<a href="#">estimate linewidths</a>	pe	get linewidth of selected peaks as width at half max
<a href="#">group peaks</a>	pg	create a peak group from selected peaks
update heights	ph	reread peak heights from NMR data file
<a href="#">label peak</a>	pl	put a label on a peak using a dialog
<a href="#">label peaks</a>	lb	put assignment labels on selected peaks
<a href="#">unlabel peaks</a>	lB	Remove labels from selected peaks
<a href="#">unoverlap labels</a>	lu	move labels so they don't overlap each other
<a href="#">right offset labels</a>	lr	move labels so they are to the right of peak markers
<a href="#">lock peaks</a>	pk	prevent peaks from moving when integrated
<a href="#">select all peaks</a>	pa	select all peaks in a spectrum
<a href="#">select fully assigned</a>	pF	select fully assigned peaks
<a href="#">select partially assigned</a>	pP	select partially assigned peaks
<a href="#">select unassigned</a>	pN	select unassigned peaks
<a href="#">select by color</a>	pC	select peaks by color
<a href="#">invert peak selection</a>	pI	invert peak selection
<a href="#">select aliased</a>	pA	select aliased peaks
<a href="#">select intra-residue</a>	pR	select i to i assigned peaks
<a href="#">select sequential</a>	pS	select i to i+1 assigned peaks
<a href="#">select medium range</a>	pM	select i to i+2, i+3 or i+4 assigned peaks
<a href="#">select long range</a>	pL	select i to i+5, i+6, ... assigned peaks
<a href="#">unintegrate peaks</a>	pU	forget that a peak has been integrated
<a href="#">unlock peaks</a>	pu	allow selected peaks to move when integrated
undo	eu	undo last peak move, deletion, or integration
<a href="#">w1 resonance peaks</a>	r1	list assigned peaks with selected peak w1 resonance
<a href="#">w2 resonance peaks</a>	r2	show assignments with selected peak w2 resonance
<a href="#">w3 resonance peaks</a>	r3	show assignments with selected peak w3 resonance

## Dialogs

<a href="#">assignment...</a>	at	make assignments
<a href="#">assignment copy...</a>	co	thresholds for copying assignments
<a href="#">assignment table...</a>	tb	table of chemical shifts
<a href="#">contour levels...</a>	ct	set contour levels and colors

<a href="#">integration...</a>	it	set integration method
<a href="#">midas...</a>	mi	display assignments on a structure using Midas
<a href="#">ornament...</a>	ot	change ornament colors and sizes
<a href="#">note...</a>	nt	add notes to peaks
<a href="#">overlay views...</a>	ol	overlay contours of one spectrum on another
<a href="#">peak list...</a>	lt	show a list of peaks for a spectrum
<a href="#">peak picking...</a>	kt	set thresholds for find peak pointer mode
pointer modes...	pm	show buttons to select pointer mode
<a href="#">predefined resonances...</a>	pd	copy resonances from another spectrum
preferences...	pf	set global Sparky options
<a href="#">print view...</a>	pt	print spectrum views to a Postscript printer
<a href="#">project view list...</a>	pv	list of views, sizes, peak counts, ...
<a href="#">python...</a>	py	Python shell for running your own programs
<a href="#">region...</a>	rt	name regions and apply operations to regions
<a href="#">rename resonances...</a>	rr	rename atoms and groups
<a href="#">resonance list...</a>	rl	list group/atom chemical shifts, standard deviations
spectrum...	st	set <a href="#">molecule and condition</a> , and <a href="#">axis shifts</a>
<a href="#">synchronize views...</a>	yt	make views scroll in parallel
<a href="#">view...</a>	vt	options controlling what is displayed in a view

## Pointer Modes

select	F1	click or drag to select or move peaks, labels, lines, ...
center	F2	clicking centers view on point
add grid both	F3	place vertical and horizontal grid lines
add grid horz	F4	add a line extending across whole spectrum
add grid vert	F5	add a line extending across whole spectrum
<a href="#">add label</a>	F6	add a text label to a peak
add line	F7	add a horizontal or vertical line between end points
<a href="#">find / add peak</a>	F8	click to place a peak or drag to find peaks
<a href="#">integrate</a>	F10	click or drag to integrate peaks
zoom	F11	drag to select a new region to show
duplicate and zoom	F12	drag to zoom and duplicate a view
<a href="#">assignment copy</a>	sh-F1	copy assignments from an already assigned spectrum
<a href="#">assignment guess</a>	sh-F2	make assignment if there is a unique guess